



- Material: ABS
- Size: 230*100*230mm
- Induction Method: Infrared
- Sensing Distance: 9cm
- Container Volume: 700mL
- Battery Power: 4AA Batteries for approximately 2100 uses
- AC Power Adapter: 100-240VAC 50/60HZ, 6VDC 2A
- 4 Sanitizer Output Options (0.5, 0.6, 0.9, and 1.3mL)

Instructions:

The Mister can be used both as a tabletop unit as well as a standing unit with the height-adjustable stand (not included in the purchase of Mister). When activating the Mister, simply stand within 15-30cm of the unit for optimum coverage.

Points of Difference:

Most Mister units use harsh compounds such as alcohol, QUATs, and other compounds as their cleaning agents. In contrast, CleanSmart uses HOCl as its main cleaning agent. Here are some points about HOCl.

HOCl:

- No Harsh Fumes
- No Harmful Chemicals
- Non-Corrosive
- Leaves no residue behind on surfaces
- No respiratory irritants
- No Allergens or asthma triggers

